


# AquaVENT<sup>®</sup>

Heated Humidified Breathing Systems



# AquaVENT® HEATER HUMIDIFIER

## PRODUCT OVERVIEW

The AquaVENT® Heater Humidifier is used to warm and humidify gases delivered to patients who require respiratory support via a tracheal tube, face mask or nasal cannula.

### The system includes:

1. Invasive or non-invasive mode selection
2. Automatic temperature selection
3. Low and high temperature alarm
4. Over-temperature protection
5. Real-time temperature tracking display allows heated plate, chamber and airway temperature to be viewed
6. Digital display
7. Servo-controlled

## PRODUCT SPECIFICATIONS

- Dimensions: W135mm x D170mm x H156mm (without chamber)
- Weight: 2.9kg (without chamber)
- Electrical Rating: Supply Frequency: 50-60Hz
- Supply Voltage: UK & EU 230V ± 23V  
NA 115V ± 12V
- Supply Current: 1.0A max at 230V
- Heater Plate Heater: 150w Heater
- Wire Capacity: 60W

## PRODUCT SPECIFICATIONS

The humidifier (excluding accessories) is warranted against manufacturing defects for a period not exceeding 2-years from the date of purchase of a previously unused unit.



## SERVICING/REPAIRS/ PREVENTATIVE MAINTENANCE

It is necessary to avoid contamination and wear of the humidifier by checking the system regularly. A review of maintenance requirements is advised once in every 12-month period. Routine maintenance will ensure continuous usage.

## ORDER CODES

### AMHH2600A-UK

AquaVENT® Heater Humidifier, servo-controlled UK

### AMHH2600A-EU

AquaVENT® Heater Humidifier, servo-controlled EU




### AMHH2600A-NA

AquaVENT® Heater Humidifier, servo-controlled NA

# AquaVENT® HEATED HUMIDIFIED BREATHING SYSTEMS

Achieve optimum delivery of High Flow Oxygen Therapy (HFOT) and Continuous Positive Airway Pressure (CPAP) with the acute care respiratory range of products from Armstrong Medical. Designed to reach the prescribed therapy quickly and maintain it consistently, these ranges allow you to focus on your patient, not the equipment, which helps to improve outcomes.

## HFOT



As open platform systems, they give you choice and control over your supply chain. It makes it easier for you to ensure the right product is always available and maintain control of your budget. Our product portfolio will easily integrate into your hospital's current supply chain.

With complete transparency on the total cost of ownership of the AquaVENT® medical equipment range you can provide the best care for your patients with confidence, knowing you've made the right choice.

## ARMSTRONG MEDICAL – IT MAKES COMMERCIAL AND CLINICAL SENSE TO CONNECT WITH US

Before you upgrade or change suppliers for your respiratory therapies consider these 8 reasons why Armstrong Medical should be part of your supply chain.

1. Complementary range of medical equipment and consumables for invasive and non-invasive ventilation therapies
2. Products designed for ease of use and straightforward, quick transition between therapies
3. BioCote® antimicrobial technology in breathing circuits provides lasting surface protection from microbial colonisation
4. Clinical training direct from the manufacturer
5. Open platform technology and compatibility with multiple devices provides flexible purchasing options and gives you control of your supply chain
6. Competitive total cost of ownership
7. Efficient inventory management
8. Increase the resilience of your supply chain in times of demand pressure



Part of



**Armstrong Medical**

Wattstown Business Park,  
Newbridge Road, Coleraine, BT52 1BS,  
Northern Ireland.

T +44 (0) 28 7035 6029

info@armstrongmedical.net



[armstrongmedical.net](http://armstrongmedical.net)

For full details of our range,  
scan QR code or visit  
[armstrongmedical.net](http://armstrongmedical.net)



Distributed by: