



**SOME  
MOISTURE  
CAN BE  
DIFFICULT  
TO MANAGE ...**

**AquaVENT® VT**  
INTUITIVE MOISTURE MANAGEMENT

**A** Armstrong  
Medical




The only circuits with BioCote<sup>®</sup> antimicrobial protection as standard

Heated inspiratory limb

Fewer filter changes





# ...BUT NOT IN AquaVENT® VT HEATED VENTILATOR CIRCUITS

Heated, vapour permeable  
expiratory  
limb maintains  
a dry circuit

Moisture in ventilator circuits can be difficult to manage, taking time and effort away from the patient.

AquaVENT® VT intuitively manages moisture throughout the system with a vapour permeable expiratory limb that allows the system to remain dry for the duration of the therapy. This allows you to focus on the patient helping to improve outcomes.

AquaVENT® VT is part of the AquaVENT® range which includes medical equipment and circuits. The entire range has been designed to allow simple and easy transition between therapies saving both time and waste. Each circuit is designed for 14 day use on a single patient.

Find full details of the AquaVENT® VT range at [www.armstrongmedical.net](http://www.armstrongmedical.net)



[armstrongmedical.net](http://armstrongmedical.net)

**Armstrong Medical**

Wattstown Business Park, Newbridge Road,  
Coleraine, BT52 1BS, Northern Ireland.

T +44 (0) 28 7035 6029

F +44 (0) 28 7035 6875

[info@armstrongmedical.net](mailto:info@armstrongmedical.net)