



OPTIMAL POSITIONING OF AquaVENT® VT VENTILATOR CIRCUIT

- Set up the circuit so that the temperature-ported connector (located at the Y-piece) is elevated above the height of the tracheal connection.
- Where an ORANGE-coloured tubing clasp is pre-assembled to the circuit, support the tubing limbs at this location.

Code	Description	Case Quantity
AMVC1775-128	<ul style="list-style-type: none"> • Adult AquaVENT® VT ventilator circuit • Dual heated limbs for use with commonly-used servo-controlled heater humidifiers • Length 1.5m with heated inspiratory limb • Heated breathable expiratory limb • 0.45m humidification limb • Bacterial viral breathing filter • Auto-fill humidification chamber 	20

