

NeoFlow[®] VT

Heated and Humidified Ventilator
Circuits for Babies



It's the little things...



NeoFlow® VT ventilator circuits are designed to maintain therapy consistently for babies.

Why NeoFlow® VT:

- Manages humidity efficiently with humidification in the inspiratory limb and vapour transmission (VT) technology in the expiratory limb, helping to prevent rainout.
- Optimum movement with the unique triple swivel Y-Piece, three-point rotation enabling smooth movement while maintaining a fixed endotracheal tube.

NeoFlow® VT ventilator circuits allow you to focus on what matters most, the baby and not the circuit.

...that make a big difference



Triple Swivel Y-Piece

Supporting care without compromise



- 1 Maintain a fixed endotracheal tube when gently repositioning baby for necessary additional care



- 2 Maintain consistent therapy with limbs that rotate independently to remove kinks



- 3 Benefit from greater mobility to safely offer bonding time with skin-to-skin care

What does Vapour Transmission (VT) technology do?

NeoFlow® VT prevents rainout and minimises condensation through a unique structure and specialised production method, allowing water vapour to escape through the expiratory limb walls, while retaining gas.




NeoFlow® VT ventilator circuits are designed to maintain therapy consistently for babies.

Why NeoFlow® VT:


- Manages humidity efficiently with humidification in the inspiratory limb and vapour transmission (VT) technology in the expiratory limb, helping to prevent rainout.
- Optimum movement with the unique triple swivel Y-Piece, three-point rotation enabling smooth movement while maintaining a fixed endotracheal tube.

NeoFlow® VT ventilator circuits allow you to focus on what matters most, the baby and not the circuit.



Rainout Prevention:

The expiratory limb lets moisture escape before condensation occurs, meaning fewer therapy disruptions.



Antimicrobial protection:

Built-in antimicrobial protection on NeoFlow® circuitry.




7-day single patient use:

Reducing waste and circuit changes.

New & Improved inspiratory limb:

New helical wire design optimises humidification and minimises rain out.



NeoFlow[®] VT

Find full details of the NeoFlow[®] VT range and product codes here



**Eakin
Healthcare**

Working together. Improving lives.



eakinrespiratory.com

Manufactured by
Armstrong Medical Limited
Wattstown Business Park,
Newbridge Road,
Coleraine, BT52 1BS,
Northern Ireland.
T +44 (0) 28 7035 6029
coleraine@eakinhealthcare.com