

AquaVENT[®]

Acute Respiratory Care
From One Trusted Supplier



The AquaVENT® range gives you everything you need to provide your chosen respiratory support for all of your patients. The entire range is designed to help you to get to the prescribed therapy quickly, and maintain it consistently.

AquaVENT® **HFOT & CPAP**
Non-Invasive Breathing Support

V2 NIV Face Masks
Patient NIV Face Mask Range

AquaVENT® **FD140i**
Dual Therapy Gas
Flow Driver

Aqua**NASE**®
HFOT Nasal Cannula

AquaVENT® **VT**
Heated Humidified
Ventilation

AquaVENT® **Bilevel**
Single Limb NIV Circuit



Helping you focus
on the patient, not
the circuit



Ease of transition
between therapies



Supported by
tailored training
from our Clinical
Engagement Team



Antimicrobial
protection as
standard on
breathing systems

Heated Humidified Ventilation

AquaVENT® VT

Heated Humidified Ventilator Circuits with a vapour permeable expiratory limb helping to prevent moisture build up in the circuit. Optimal humidification, dry expiratory limb, BioCote® anti-microbial protection.

HFOT & CPAP

AquaVENT® HFOT & CPAP Breathing Support

Circuits designed to transition easily between therapies, supplying only what you need, reducing waste and saving money. Easy transition prioritises patient wellbeing. BioCote® anti-microbial protection.

HFOT Nasal Cannula

AquaNASE® Nasal Cannula

Nasal interface for effective delivery of HFOT with a double head strap for enhanced stability and comfort. The lanyard securement supports the weight of the breathing circuit to prevent prong displacement and therapy disruption.

Acute Respiratory Care

Everything you need to provide the best care for your patients in the acute setting, from one trusted supplier.

AquaVENT® FD140i

Designed with clinicians and patients in mind, the AquaVENT® FD140i dual therapy gas flow driver delivers safe and effective HFOT and CPAP therapies. It provides an alternative option to using an invasive ventilator to deliver non-invasive therapies when facilitating patient weaning, offering simple transitions between therapies.

Dual Therapy Gas Flow Driver

V2 NIV Face Masks

The V2 NIV face mask range has a unique design to maximise patient comfort and minimise leaks to ensure optimal therapy delivery. The range includes options with CO₂ venting, non-vented, universal (with detachable CO₂ venting) and with or without anti-asphyxia valve.

Patient NIV Face Masks

AquaVENT® Bilevel

Our versatile range of non-invasive, single-limb NIV circuits available with an optional supplemental oxygen line, a bacterial/viral filter, and a CO₂ exhalation valve, are suitable for use with most bilevel or mechanical ventilators, offering clinicians flexible options to meet diverse patient needs.

Single Limb NIV Circuits

Supporting Respiratory Products

Extensive range of respiratory support products



**CombiFlex®
Circuits**



Ultra-PEP®

Mapleson C



Benefit from one trusted supplier for all your patients' respiratory care needs. With our wide product range, we offer solutions tailored to your practice - delivering optimal therapy quickly, effectively, and comfortably for every patient.

Explore our full range on our website and contact our team for personalised support.



eakinrespiratory.com

Manufactured by:
Armstrong Medical Limited
Wattstown Business Park,
Newbridge Road,
Coleraine, BT52 1BS,
Northern Ireland.
T +44 (0) 28 7035 6029
coleraine@eakinhealthcare.com